



PVDOMICS STUDY

Event Information Sent to the DCC Form #402

Instructions: The Data Coordinating Center (DCC) will notify clinical center staff to complete and enter Form 402 when an event (hospitalization, ER visit, other SAE and/or death) shows that a packet needs to be scanned and uploaded to DCC folder in PVDOMICS via SSH. See the MOP for detailed instructions on processing the packet.

Form 401 must be entered by the clinical center before this form is entered into the database.

NOTE: Do NOT send any packets to the DCC unless notified to do so by the DCC.

--	--	--	--	--	--

1. Identification Number

--	--

2. Alphacode

--	--	--	--	--	--	--	--

3. Date of onset of event (mm/dd/yyyy)

4. Date event packet scanned and transferred to the DCC? (mm/dd/yyyy)..... ____ / ____ / ____

5. Type of information scanned and transferred to the DCC:

- a. Discharge summary (0=No, 1=Yes) ____
- b. ER summary note (0=No, 1=Yes) ____
- c. Physician's narrative summary (0=No, 1=Yes)..... ____
- d. Expiration summary (0=No, 1=Yes)..... ____
- e. Autopsy report (0=No, 1=Yes) ____
- f. Death certificate (0=No, 1=Yes) ____
- g. Other documentation sent (0=No, 1=Yes) ____

If other, describe other documentation provided:

200. Date form completed (mm/dd/yyyy)..... ____ / ____ / ____

201. Username of PVDOMICS Study staff member completing this form..... ____

Clinical Center Use Only

Date form entered (mm/dd/yyyy) ____ / ____ / ____

Username of person entering this form ____